

培养基配方表 (Media Formulation)

名称: Hank's 平衡盐溶液 HBSS 含碳酸氢钠和钙、镁离子, 不含酚红

货号: B410KJ

组分 (Components)	分子量 (MW)	浓度 (mg/L)	摩尔浓度 (mM)
Potassium Chloride	74.55	400	5.3333335
Potassium Phosphate monobasic	136.1	60	0.44117647
Sodium Bicarbonate	84.01	350	4.1666665
Sodium Chloride	58.44	8000	137.93103
Sodium Phosphate dibasic	142	48	0.33802816
Calcium Chloride	111	140	1.2612612
Magnesium Chloride hexahydrate	203.3	100	0.49261084
Magnesium Sulfate	120.4	48.8	0.40650406
D-(+)-Glucose	180.2	1000	5.5555553